

MARIA® for Foods - Validated Performance Parameters

| Allergen: | Ana o 3 | Ara h 1 | Ara h 3 | Ara h 6 | Api g 1 | Bos d 5 (Native) | Bos d 11 | Cor a 9 | Cyp c 1 |
|--|-----------|-----------|----------|-----------|-----------|---------------------|-----------|-----------|----------|
| Linearity (R²)¹ | 0.998 | 0.996 | 0.999 | 0.997 | 0.995 | 0.998 | 0.997 | 0.999 | 0.995 |
| Range (ng/ml)² | 20-0.01 | 250-0.12 | 125-0.06 | 12.5-0.02 | 200-0.10 | 50-0.10 | 500-0.98 | 50-0.02 | 250-0.12 |
| Limit of Quantification³ | | | | | | | | | |
| <i>LLOQ (ng/ml)^{3a}</i> | 0.01-0.02 | 0.12-0.98 | 0.06 | 0.02 | 0.20-0.78 | 0.10 | 0.98-3.91 | 0.02-0.10 | 0.12 |
| <i>ULOQ (ng/ml)^{3b}</i> | 20-5 | 250 | 125 | 50-6.25 | 200-100 | 50 | 500-250 | 50 | 250-62.5 |
| Accuracy (% Recovery)⁴ | | | | | | | | | |
| <i>Intra-assay (n=9)^{4a}</i> | 101-119% | 92-113% | 89-101% | 91-102% | 83-106% | 82-115% | 89-132% | 78-94% | 89-132% |
| <i>Inter-assay (n=36)^{4b}</i> | 111% | 101% | 94% | 97% | 94% | 92% | 106% | 85% | 113% |
| Precision (%CV)⁵ | | | | | | | | | |
| <i>Intra-assay (n=9)^{5a}</i> | 2-6% | 1-13% | 4-9% | 3-6% | 2-10% | 4-6% | 3-15% | 3-6% | 3-15% |
| <i>Inter-assay (n=36)^{5b}</i> | 5% | 6% | 7% | 5% | 7% | 6% | 7% | 5% | 6% |

| Allergen: | Gal d 1 | Gal d 2 | Gly m 5 | Jug r 1 | Pru du 6 | Ses i 1 | Shrimp Tropomyosin | Sin a 1 |
|--|-----------|---------|-----------|---------|----------|-----------|-----------------------|-----------|
| Linearity (R²)¹ | 0.992 | 0.997 | 0.998 | 0.999 | 1.000 | 0.998 | 0.997 | 0.999 |
| Range (ng/ml)² | 500-0.98 | 50-0.02 | 500-0.98 | 20-0.01 | 125-0.24 | 12.5-0.02 | 25-0.02 | 15-0.01 |
| Limit of Quantification³ | | | | | | | | |
| <i>LLOQ (ng/ml)^{3a}</i> | 0.24-1.95 | 0.02 | 0.24-0.49 | 0.01 | 0.24 | 0.02-0.05 | 0.02 | 0.01-0.12 |
| <i>ULOQ (ng/ml)^{3b}</i> | 500-250 | 50-12.5 | 500-125 | 20 | 125 | 12.5-6.25 | 50-12.5 | 15 |
| Accuracy (% Recovery)⁴ | | | | | | | | |
| <i>Intra-assay (n=9)^{4a}</i> | 83-97% | 97-118% | 95-111% | 80-131% | 71-102% | 79-117% | 94-113% | 72-131% |
| <i>Inter-assay (n=36)^{4b}</i> | 91% | 106% | 106% | 98% | 90% | 96% | 103% | 97% |
| Precision (%CV)⁵ | | | | | | | | |
| <i>Intra-assay (n=9)^{5a}</i> | 4-8% | 2-7% | 3-9% | 0-12% | 4-15% | 2-14% | 4-6% | 2-22% |
| <i>Inter-assay (n=36)^{5b}</i> | 6% | 4% | 6% | 5% | 7% | 6% | 5% | 8% |

Validation Notes:

1. **Linearity** is the mean R^2 of four MARIA plates for control curves generated using 4-parameter logistic fit.
2. **Range** is the average usable range of control curves from four MARIA plates that have a value of 70-130% of the expected concentration, with %CV < 15 between duplicate points.
3. **Limit of Quantification**
 - 3a. LLOQ - The lowest concentration points of four control curves with a recovery of 70-130% and %CV <15, expressed as a range.
 - 3b. ULOQ - The highest concentration points of four control curves with a recovery of 70-130% and %CV <15, expressed as a range.
4. **Accuracy**
 - 4a. Intra-assay - The range of average percent recovery of samples A, B, and C run in triplicate from four MARIA plates (n=9).
 - 4b. Inter-assay - The overall average percent recovery of samples A, B, and C run in triplicate from four MARIA plates (n=36).
5. **Precision**
 - 5a. Intra-assay - The range of average percent coefficient of variation of samples A, B, and C run in triplicate from four MARIA plates (n=9).
 - 5b. Inter-assay - The overall average percent coefficient of variation of samples A, B, and C run in triplicate from four MARIA plates (n=36).