

Quantitative ELISA kits and calibrated allergen standards

A comprehensive range, cited in over 400 scientific publications

Each ELISA kit contains:

- capture monoclonal antibody
- calibrated allergen standard
- calibrated biotinylated secondary mAb (or rabbit polyclonal antibody)*
- sufficient reagents for ten 96 well micro-titer plates (ELISA-10)
- full assay protocol and references supplied

Quantity discounts on ELISA kits

- 5% discount for 3 kits
- 10% for between 4 and 9 kits
- 15% for more than 10 kits

Quantitative ELISA kits

Code	Description	Antibody Pairs	Quantity	
House dust mite				
EL-DP1	Der p 1 ELISA-10	5H8/4C1	1 Kit	
EL-DP1A	Der p 1 ELISA-10	10B9/5H8	1 Kit	
EL-DF1	Der f 1 ELISA-10	6A8/4C1	1 Kit	
EL-DF2	Der f 2 ELISA-10	1D8/7A1	1 Kit	
EL-D2	Mite Group 2 ELISA-10	1D8/7A1	1 Kit	
EL-BT5	Blo t 5 ELISA-10	4D9/4G9	1 Kit	
Animal				
EL-FD1	Fel d 1 ELISA-10	6F9/3E4	1 Kit	
EL-CF1	Can f 1 ELISA-10	10D4/R anti Can f 1	1 Kit	
EL-RN1	Rat n 1 ELISA-10	RUP-6/RUP-1	1 Kit	
EL-MM1	Mus m 1 ELISA-10	R anti Mus m 1	1 Kit	
EL-EC4	Equ c 4 ELISA-10	1O3/14G4	1 Kit	
EL-BD2	Bos d 2 ELISA-10	3D4/mAb1	1 kit	
Cockroach				
EL-BG1	Bla g 1 ELISA-10	10A6/R anti Bla g 1	1 Kit	
EL-BG2	Bla g 2 ELISA-10	7C11/R anti Bla g 2	1 Kit	

*Not supplied: labelled goat anti-rabbit antibodies; streptavidin peroxidase; ABTS color developer; routine laboratory chemicals. See Protocols on www.inbio.com for full details.

Code	Description	Antibody Pairs	Quantity	
Mold				
EL-AF1	Asp f 1 ELISA-10	4A6/R anti Asp f 1	1 Kit	
EL-AA1	Alt a 1 ELISA-10	2C10/3B6	1 Kit	
EL-SCX	<i>S. chartarum</i> X ELISA-10	3E11/2C9	1 Kit	
EL-SCY	<i>S. chartarum</i> Y ELISA-10	3F10/3F10	1 Kit	
EL-AVX	<i>A. versicolor</i> ELISA-10	5F5/9A5	1 Kit	
Pollen				
EL-BV1	Bet v 1 ELISA-10	4B10/2E10	1 Kit	
EL-PP5	Phl p 5 ELISA-10	1D11/B01	1 Kit	
EL-AM1	Amb a 1 ELISA-10	R anti Amb a 1/R anti Amb a 1	1 Kit	
Food				
EL-AH1	Ara h 1 ELISA-10	2C12/2F7	1 kit	
EL-AH2	Ara h 2 ELISA-10	1C4/R anti Ara h2	1 kit	
EL-AH6	Ara h 6 ELISA-10	3B8/3E12	1 kit	
* Latex - FITkit®				
EL-HB1	FITKit® Hev b 1 ELISA-1		1 kit	
EL-HB3	FITKit® Hev b 3 ELISA-1		1 kit	
EL-HB5	FITKit® Hev b 5 ELISA-1		1 kit	
EL-HB6	FITKit® Hev b 6.02 ELISA-1		1 kit	

*All inclusive kit: FITKits® contain 1 pre-coated 96 well microtiter plate and all necessary reagents. Indoor Biotechnologies is a distributor of FITkit® for latex allergens, produced by Icosagen AS, Tartu, Estonia.

Calibrated allergen standards for ELISA

Code	Description	Quantity	
Universal allergen standard			
ST-UAS	Universal allergen standard	1 vial	
House dust mite			
ST-BT5	Blo t 5	1000ng	
ST-DF2	Der f 2	400µl	
Cockroach			
ST-BG1	Bla g 1	10 units	
Mold			
ST-AF1	Asp f 1	160ng	
ST-AA1	Alt a 1	400ng	
ST-AV1	Ave x 1	400ng	
ST-SC1	Sch X and Sch Y	400ng	