

Product Code: NA-CA14-1

Gel Image:
nCor a 14

Allergen: nCor a 14 (*Corylus avellana* allergen 14)

Lot No: XXXXX

Source: 2S Albumin from Hazelnut flour.

Mol. Wt: 10 kD

Purification: Purified from hazelnut extract by Immunoaffinity & HPLC.
Purity >95% by SDS-PAGE.

Concentration: See Product Insert.

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See Product Insert.



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com