

Product Code: NA-CF3-1

Allergen: nCan f 3 (*Canis familiaris* allergen 3)

Lot No: XXXXX

Source: Partially purified dog albumin

Mol. Wt: 69 kD

Purification: Purified by multi-step column chromatography.
Purity >95% by SDS-PAGE.

Concentration: See product insert.

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Notes: (1) Can f 3 is a serum albumin

Gel Image:
nCan f 3



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com