

Product Code: NA-CJ1-1

Allergen: nCry j 1 (*Cryptomeria japonica* allergen 1)

Lot No: XXXX

Source: Japanese cedar pollen (*Cryptomeria japonica*)

Mol. Wt: 40 kD

Purification: From Japanese cedar pollen extract by multi-step chromatography.
Purity > 90 % by SDS-PAGE

Concentration: 1.1 mg/ml by OD 280 nm.
Extinction Coefficient = 1.6

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Notes: (1) Purified natural Cry j 1 runs at 40 kD with dimer and trimer on non-reducing SDS-PAGE.

Gel Image:

nCry j 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com