

Product Code: NA-GM6-1

Gel Image:

Allergen: nGly m 6 (Glycine max allergen 6)

nGly m 6



Lot No: XXXXX

Source: Glycinin from soybean.

Mol. Wt: Multiple subunits ranging from 37-50 kD.

Purification: Purified from soybean flour by multi-step HPLC.
Purity > 95 % by SDS-PAGE.

Concentration: See product insert.

Formulation: Preservative and carrier-free in 1X PBS pH 7.4
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com