

Product Code: RE-AH8-1

Allergen: rAra h 8 (*Arachis hypogaea* allergen 8.0101)

Lot No: XXXX

Source: *E. coli*

Mol. Wt: 17.8 kD

Purification: Purified from *E. coli* culture by HPLC gel filtration.
Purity > 95 % by SDS-PAGE.

Concentration: See product insert.

Formulation: 20mM Tris, 0.3M NaCl, pH 8.0. Filtered through 0.22µm filter

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

Notes: (1) rAra h 8 has a C-terminal 6xHis-tag.

Gel Image:

rAra h 8



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com