

Product Code: RE-CS5-1

Allergen: rCan s 5 (*Cannabis sativa* allergen 5.0101, Indian hemp)

Lot No: XXXX

Source: *E. coli*

Mol. Wt: 18 kD

Purification: Multi-step chromatographic purification.
Purity >95% by SDS-PAGE.

Concentration: 1.3 mg/ml by OD 280 nm; Extinction coefficient = 0.647.

Formulation: Preservative-free and carrier-free in 20mM Tris, 0.3M NaCl, pH 8.0.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXX; when stored at -20°C.

Notes: (1) rCan s 5 has a N-terminal 6x His-tag.

Gel Image:
rCan s 5



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com