

Product Code: RP-AA1-1

Allergen: rAlt a 1 (*Alternaria alternata* allergen 1)

Lot No: XXXX

Source: *Pichia pastoris*

Mol. Wt: ~15kD

Purification: Purified from *Pichia pastoris* culture by HPLC gel filtration. Purity >95% by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.


Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rAlt a 1 produced in *Pichia pastoris* has been observed as a dimer at ~30kD in the SDS-PAGE gel.

Gel Image:

rAlt a 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com