

Product Code: RP-BG1-1

Allergen: rBla g 1.0101 (*Blattella germanica* allergen 1)

Lot No: XXXX

Source: *Pichia pastoris*

Mol. Wt: Multiple molecular forms.
Purity >95 % by SDS-PAGE.

Purification: Purified from *Pichia pastoris* culture by affinity chromatography.

Concentration: See product insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rBla g 1 is produced in *P. pastoris*. As a result of proteolytic cleavage of tandem amino acid repeats, the product contains multiple molecular forms of rBla g 1.

Gel Image:
rBla g 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com