

Product Code: RE-BG5-1

Allergen: rBla g 5.0101 (*Blattella germanica* allergen 5)

Lot No: XXXX

Source: *E. coli*

Mol. Wt: 23 kD

Purification: Purified from lysed bacterial supernatant over glutathione Sepharose. Purity > 95% by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.


Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) Bla g 5 belongs to the glutathione transferase protein family.

Gel Image:

rBla g 5



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com