

Product Code: RP-CA1-1

Allergen: rCor a 1 (*Corylus avellana* allergen 1.0401)

Lot No: XXXX

Source: *Pichia pastoris*

Mol. Wt: 18 kD

Purification: Purified from *Pichia pastoris* culture by His-Trap chromatography. Purity > 95 % by SDS-PAGE, run under reducing conditions.

Concentration: See product insert

Formulation: Preservative and carrier-free in 1X PBS, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rCor a 1 has a C-terminal 6xHis-tag.
(2) Substituted Q for N (N126Q) at Glycosylation site 126.

Gel Image:

rCor a 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com