

Product Code: RP-CA8-1

Allergen: rCor a 8 (*Corylus avellana* allergen 8.0101)

Lot No: XXXX

Source: nLTP from hazelnut expressed in *Pichia pastoris*

Mol. Wt: 9 kD; under reducing conditions

Purification: Purified from *Pichia pastoris* culture by multi-step HPLC.
Purity > 95 % by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative and carrier-free in PBS, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rCor a 8 runs as doublet.
(2) rCor a 8 has a C-terminal 6xHis-tag.
(3) A putative N-glycosylation site was silenced (N33Q).

Gel Image:

rCor a 8



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com