

Product Code: RP-CF2-1

Allergen: rCan f 2 (*Canis familiaris* allergen 2.0101)

Lot No: XXXX

Source: *Pichia pastoris*

Mol. Wt: 19 kD

Purification: Purified from *Pichia pastoris* culture by multi-step HPLC.
Purity > 95 % by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.


Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rCan f 2 presents as a 19kD unglycosylated and a 23 kD glycosylated monomer. In addition, a dimer around 40kD is visible under non-reducing conditions.
(2) rCan f 2 has a C-terminal 6xHis-tag

Gel Image:

rCan f 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com