

Product Code: RP-PP3-1

Allergen: rPru p 3 (*Prunus persica* allergen 3.0102)

Lot No: XXXX

Source: nLTP from peach expressed in *Pichia pastoris*

Mol. Wt: 10 kD

Purification: Purified from *Pichia pastoris* by multi-step HPLC.
Purity > 95 % by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative and carrier-free in PBS, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rPru p 3 has a C-terminal 6xHis-tag.

Gel Image:

rPru p 3



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com