

Product Code: RP-RN1-1

Allergen: Recombinant Rat n 1.0101

Lot No: XXXX

Source: *Pichia pastoris*

Mol. Wt: 19.5 kD

Purification: Purified by Nickel IMAC chromatography.
Purity > 95% by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) To avoid hyper-glycosylation, N-glycosylation motif at position 35
was mutated (N35Q).
(2) rRat n 1.0101 has a C-terminal 6xHis-tag.

Gel Image:

rRat n 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com