

Product Code: RPI-AV3-1

Allergen: rArt v 3 (*Artemisia vulgaris* allergen 3.0101)

Lot No: XXXX

Source: LTP from mugwort expressed in *Pichia pastoris*

Mol. Wt: 11.6 kD

Purification: Multi-step chromatographic purification.
Purity >95% by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rArt v 3 has a N-terminal 8x His-tag.

Gel Image:
rArt v 3



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com