

Product Code: RPI-JR3-1

Allergen: rJug r 3 (*Juglans regia* Allergen 3.0101)

Lot No: XXXX

Source: *Pichia pastoris*

Mol. Wt: 10 kD

Gel Image:

rJug r
3



Purification: Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% on silver stained SDS-PAGE.

Concentration: See Product Insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, 300mM NaCl, pH 7.4. Sterile filtered.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See Product Insert

Notes: (1) rJug r 3 has a N-terminal 6xHis-tag.
(2) Made in USA.

For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com