

**Product Code:** RPI-PA1-1

**Allergen:** rPen a 1 (*Penaeus aztecus* allergen 1.0101)

**Lot No:** XXXXX

**Source:** *Pichia pastoris*

**Mol. Wt:** 33 kD

**Purification:** Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% by SDS-PAGE.

**Concentration:** See Product Insert.

**Formulation:** Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.


**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See Product Insert; when stored at -20°C.

**Notes:** (1) rPen a 1 has a C-terminal 6xHis-tag.

**Gel Image:**

rPen a 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)