

Product Code: RPI-PL1-1

Allergen: rPla I 1 (*Plantago lanceolata* allergen 1.0101)

Lot No: XXXXX

Source: *Pichia pastoris*

Mol. Wt: 15 kD

Purification: Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% by SDS-PAGE.

Concentration: See Product Insert.

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See Product Insert; when stored at -20°C.

Notes: (1) rPla I 1 has a N-terminal 6xHis-tag.

Gel Image:

rPla I 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com