

LoTox™ Egg-White Protein

Product Code: LTN-GDE-1


Allergen: *Gallus domesticus (G. gallus)*

Lot No: XXXXX

Source: Natural egg-white powder

Mol. Wt: Gal d 1: 28kD, Gal d 2: 44kD, Gal d 3: 78kD, Gal d 4: 14kD

Gel Image



Egg-white Protein

LC-MS/MS: See product insert.

ELISA: See product insert.

Concentration: See product insert.

Formulation: Preservative and carrier-free in sterile, 0.5M NaCl, 0.1M Potassium Phosphate pH 7.40. Sterile filtered.

Storage: Store at -20°C. Re-freeze unused portion and aliquot into smaller quantities to avoid freeze/thaw cycles. Gently defrost at room temperature.

Expiration: See product insert.

Released:

Lixin D Vaile

Notes:

- (1) (1) A LoTox™ product: Endotoxin ≤ 0.03 EU/ μ g.
- (2) Relative abundance results are based on label-free quantification of precursor peptide ion intensity (peak area, incl. unique + razor peptides) using Proteome Discoverer 2.2 algorithms.
- (3) Calculations of relative abundance can be variable and depend on multiple factors such as sample complexity, instrument resolution and analysis software.

For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA 22903
United States
(434) 984-2304
mail@inbio.com

Vision Court, Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom
+44 (0) 2921 674 640
Info@inbio.com