

October, 2022

Re: Bovine Serum Albumin (BSA)


Effective immediately, InBio recommends the use of Bovine Serum Albumin (BSA) Fraction V manufactured by Roche Diagnostics (catalog number 03116964001) for ELISA 2.0 and MARIA® assay buffers. Buffers prepared with BSA from other manufacturers may become cloudy even after sterile filtration, leading to increased MARIA® analysis time and a higher incidence of bead aggregation errors.

In accordance with Good Laboratory Practices, InBio recommends monitoring the quality of prepared buffers and advises that assay results may be affected if using BSA supplied by manufacturers other than Roche Diagnostics.

Please note that InBio has no financial conflict or other interest in Roche Diagnostics. Our recommendation is based solely on product performance.

Please contact Bryan Smith, Research Scientist and Quality Assurance Specialist, at InBio with any questions or comments regarding this matter.

Bryan Smith



Research Scientist and Quality Assurance Specialist
InBio
bsmith@inbio.com