

## Food

## Inhaled

Target Allergen	hIgE mAb
<b>Peanut</b>	
Ara h 1	3B10 4G4
Ara h 2	16A8 9H11 13D9 11F10 26C3
Ara h 3	3C3 14G12*
Ara h 6	8F3 20G11 7B6 15C2

Target Allergen	hIgE mAb
<b>Red Meat</b>	
Alpha-gal	16D9 10H8
<b>Milk</b>	
Bos d 8	1D5
<b>Egg</b>	
Gal d 2	11B6
Gal d 4	1E7
<b>Walnut</b>	
Jug r 1	9H7
<b>Cashew</b>	
Ana o 3	2F5

Target Allergen	hIgE mAb
<b>Dust mite</b>	
Der p 2	2G1 2F10 1B8
<b>Cat</b>	
Fel d 1	11A12 6A1 1B7
Fel d 4	1E18
<b>Dog</b>	
Can f 1	1J11

\*Not reactive with rAra h 3 ImmunoCAP (f424), reactivity confirmed by biotinylated nAra h 3 streptavidin cap.